

CARESS Document Set

User Commands

Busse/Brehmer/Sauer

Printed on December 13, 2013
SVN-Revision: 1701 (2010-09-01)

Abstract: This document, part 2 of the CARESS document set, lists all the CARESS user commands. It is both the command reference and the command user's guide and has been updated to reflect the changes in the current CARESS release.

Users new to CARESS will find here descriptions of the command's structure and purposes. All available commands are listed with some notes about the meaning and a 'real world' example.

Contents

1	Overview	1
2	Command Structure	2
3	Commands	5
4	General setup commands	6
4.1	Definition of an User Experiment	6
4.2	Definition of Names, Title, Comment, etc.	6
4.3	Parameters for Driving Modules	7
4.4	Definition of Acquisition Parameters	7
4.5	Detector related parameters	8
4.6	Define Online Protocol and Plot of Scan Data	9
4.7	Reset internal setup lists	10
5	Special setups	12
5.1	Polarisation	12
5.2	Coordinate Transformations (Angles \Leftrightarrow H,K,L, ...)	13
5.3	Diverse special	14
6	Data acquisition and other actions	16
6.1	Perform Basic Actions	16
6.2	Scans	16
6.3	Initialisation	18
7	Special actions	19
7.1	Reference movements	19
7.2	Universal device access	19
7.3	Complex actions	21
8	Control special hardware devices	22
8.1	Definition of Compound Positions, etc.	22
8.2	Define Data Acquisition Unit DAU	22
8.3	Neutron Chopper System	23
9	Flow control	24
9.1	Wait commands	24
9.2	Input Utilities	25

9.3	Control CARESS by an external user program, ...	26
9.4	Interrupt / Kill / Continue / Terminate CARESS Activities	27
10	Miscellaneous commands	29
Index		31

1 Overview

Commands serve different purposes:

General setup define an user experiment, names, titles, comments and parameters

Special setups dealing with polarisation, coordinate transformations and some others

Data acquisition performing basic actions und run scans

Special actions like reference movements and universal device access

Special hardware access to deal with complex hardware systems

Flow control do loops and control the behaviour

Miscellaneous diverse useful pieces

2 Command Structure

CARESS command input is requested by the prompt

```
CAR>
```

in case of exceptional situations with

```
EXC>
```

An input line may contain only one command, but a command may spread several lines. If the last character is '-', a continuation line is expected. Comment lines start with ';'. A comment starts with the command *COM*. A comment on a command line starts with '!'.
Any sequence of CARESS commands can be read from a text file generated using an editor or the CARESS GUI. Input from such a command file (e.g. `mycommands.com`¹) is invoked by

```
CAR>@mycommands
```

Command files may be nested.

A command starts with a command key which is optionally followed by a parameter list, the parameter list may be

- a list of values, e.g. `DS 3.3542 3.3542`
- a list of subkeys (mostly module names) with or without values, e.g.

```
TOL OMGS=0.02,PHIS=0.01  
SL1 TIM1,SDET
```

Some² commands are *additive*, i.e. information in several occurrences of the same command is added (merged) e.g.

```
BLS OMGS=-100 100  
BLS PHIS=0 75
```

have the same effect as

¹Note: the suffix *.com* is mandatory.

²valid for: *BLS, DEFOFF, DEFON, DEFTURN, PRW, PRWA, SOF, TMO, TOL, WIND, WINDA, SELDET, DEFPOS, DEFINE, DEFSAM*

```
BLS OMGS=-100 100 PHIS=0 75
```

All other commands are not additive, and the information of the last occurrence of such a command is relevant, only.

To obtain the current coverage of a command (with no matter if additive or not) it is usefull to call the CARESS show command.

For the example above:

```
CAR> SHOW BLS
```

yields to this message on the protocol window:

```
BLS OMGS=-100 100 PHIS=0 75
```


3 Commands

This chapter lists all the CARESS commands. Which commands are accepted at a specific instrument (or in a specific situation) depends on the functionality and setup.

Typographical notes concerning the following command definitions:

{...}	means <i>repetition, at least one occurrence</i>
[...]	means <i>optional</i>
<real>	::= PASCAL notation, in addition: {{<digit>}}{<digit>}
<identifier>	::= <letter>{<letter or digit or _>}, max. 20 characters
<text>	::= {<symbol>}, max. 68 characters
<symbol>	::= <program name>
<string>	::= <identifier> '<text>'
<file_name>	::= <identifier>
<peakinfo>	::= ON OFF LL <beta><delta_b><peak form factor> NLL <left points><right points>

Each command is listed with:

- command purpose (short)
command syntax
command example

4 General setup commands

Most of the commands of this group are essential for all of the supported instruments.

4.1 Definition of an User Experiment

- Instrument selection information

EXPTYPE *<instrument>**<axes or mode>**<detector>**<exptype>**<supplementary equipment>*

EXPTYPE E1 A3 SDET EULER TEMP

EXPTYPE V1, A2, ADET, POWD

- List of modules depending on EXPTYPE command

ONL {*<module>*}

ONL OMGS TTHS TEMP

4.2 Definition of Names, Title, Comment, etc.

- User Name

USN *<identifier>*

USN MILLER

- Title

TITLE *<text>*

TITLE This is a title.

- Sample Name

SAM *<identifier>*

SAM CU

- List of samples (for V4)

DEFSAM {*<pos_no>**<string>**<value>**<value>*}

if sample changer

DEFSAM 1 'HMI#005' 0.0 0.0, 2 B60 0.0 0.0

```
otherwise
DEFSAM 1 'HMI#005' 0.0 0.0
```

- Comment

```
COM <text>
```

```
COM This is a comment.
```

4.3 Parameters for Driving Modules

- Limits (Range)

```
BLS {<module> = <lower limit>, <upper limit>}
```

```
BLS OMGM=0,100 TTHS= -50,70 OMGS=40,80
```

- Offset (real position = encoder_value + offset_value)

```
SOF {<module> = <real>}
```

```
SOF TTHM=0 OMGS=20.5
```

- Tolerance

```
TOL {<module> = <tolerance value>}
```

```
TOL OMGS=0.02 PHIS=0.01
```

- Timeout in seconds (for devices with known velocity: supplementary timeout)

```
TMO {<module> = <timeout value>}
```

```
TMO OMGS=10, TEMP=200
```

4.4 Definition of Acquisition Parameters

- Master (preset) Counter

```
MMI <module> <preset value>
```

```
MM1 MON 10000
```

- Slave (acquisition) Modules

```
SLI {<module>}
```

```
SL1 TIM1, TTHS, LDET, TEMP
```

- Protocol Modules (read at begin and at termination of (partial) data acquisition)

```
PRO {<module>}
```

```
PRO TEMP
```

- Definition of Supplementary Modules to be journalised like a SET Module during a Scan.

```
SCANV {<module>}
```

```
SCANV MAGF
```

- Definition of Permanent Movement during Acquisition

```
ACQPMOVE <module> = <pos1><pos2>
```

```
ACQPMOVE OMGS=-100 100
```

4.5 Detector related parameters

- The user can SElect one or more (up to 4) DETectors out of the single detectors and/or multi detector for special purpose (define cell or region of cells). This command applies to the Data Acquisition Unit (DAU) only. The predefined names for selected detectors are D1, D2, D3 and D4. You can omit parameters, they are assumed as zero. For a single detector the y_chan should be zero in any case.

```
SELDET {Di = <x_chan>[<y_chan>[<radius>]]}
```

select 4 cells of a linear detector:

```
SELDET D1=17 D2=20 D3=21 D4=2
```

select a region: cells 14 ... 20 of a 1-dim multi detector:

```
SELDET D1= 17, 0, 3
```

select a region of an area detector:

```
SELDET D1 = 10, 20, 5
```

- Erase the selection of detector cell or region of cells of a detector.

```
RESET SELDET
```

```
RESET SELDET
```

- Up to 4 Windows for Linear Multi Counter (i = 1 ... 4)

```
WIND {Wi = <lower channel>, <upper channel>}
```

```
WIND W1=20 80 W2=150 250
```

- Up to 4 Windows for Area Detector (i = 1 ... 4)

```
WINDA {Wi = <lower_x_chan>, <upper_x_chan>, <lower_y_chan>, <upper_y_chan>}
```

```
WINDA W1 = 20 80 60 128
```

- Correction of Multi Detector Counts (CARESS will append extension .dat to the filename)

```
CALIB <module> = <file_name>
```

```
CALIB LDET = e2_ldet_corr
```

- Correspondence Detector Angle \Leftrightarrow channel number of multi counter

RELA {*module*} = *channels*, *ref_channel*, *cell_distance*} –
[, *distance sample – detector*] = 0]

1-dim multi detector with fixed distance sample-detector, *cell_distance* in degrees:

RELA TTHS = 400 1 0.2

area detector with variable distance sample-detector, *cell_distance* in cm:

RELA TTHS = 128 1 0.1, NYS=128,1,0.1, YSD=0

Note:

RELA sets the parameters for the detector given by *EXPTYPE* (refer to page 6). If more than one area detector is involved, use *MRELA*.

- Correspondence Detector Angle \Leftrightarrow channel number of multi counters

MRELA {*detector*}{*module*} = *ref_channel*, *offset*}

two 2-dim multi detectors, absolute distances:

MRELA ADET1 TTHS=128 0, NYS=128,0, YSD=0,0

MRELA ADET2 TTHS=128 0, NYS=128,0, YSD=0,0

two 2-dim multi detectors, relative distances:

MRELA ADET1 TTHS=128 64, NYS=128,128, YSD=0,100

MRELA ADET2 TTHS=128 TTH1, NYS=128,64 YSD=0,100

Note:

Each area detector is to define with one *MRELA* command. A complete setup includes a corresponding entry in *corbadevice.dat* (*startcommands* += *MRELA*).

The offset is relative to the reference channel of the module.

4.6 Define Online Protocol and Plot of Scan Data

- Protocol Parameters

PRT ON/OFF/FUL/SHT

FUL: full, SHT: short

PRT FUL

- Define modules to be printed in the protocol terminal

PRTMOD {*device*}

PRTMOD OMGX, RESIST, MON2

Note:

This command has no effect on the data file it just reduces the number of *online* printed protocol modules to fit screen size etc. The files contain all *SLI* and *PRO* modules.

- Define Protocol Windows one dimensional (in case of SDET or LDET)

PRW {*Wi* = *<lower_channel><upper_channel>*}
(i=1 ... 4)

PRW W1= 5 196

two dimensional (in case of ADET)

PRWA {*Wi* = *<lower_x_ch><upper_x_ch><lower_y_ch><upper_y_ch>*}
(i = 1 .. 4)

PRWA W1 = 1 128 1 80

- Peak calculation

LL: method lehman-larsen;

NLL: not LL

PEAK *<peakinfo>*

PEAK ON

PEAK LL, -1, -1, 0.5

PEAK NLL 8, 5

- Switch OFF/ON Online Plot, define Refresh Time (in seconds)

DISP OFF/ON/LIV

DISP LIV 10

- Define Histogram Axes

SPEC [*<x_axis>*]*<y_axis>*

<x_axis> : any scan module,

<y_axis> : any slave module

SPEC OMGA SDET

SPEC LDET

- Supervise environment modules by saving their values periodically

LOGDEV PERIOD = *<seconds>*{*<environment module>*}

(save action not yet implemented)

LOGDEV PERIOD=60 TEMP, MAGF

4.7 Reset internal setup lists

- Reset / Delete Command Parameters

with commands: CALIB, SELDET, PERMOVE, PASS, TOF, SCANV, WAIT1, WAIT2, ONL, OFL, SL1, PRO, WIND, WINDA, PRW, PRWA, RESOSTEP, SPEC

RESET {<commandkey>}

RESET PASS

5 Special setups

The following commands will apply to some special configurations at some instruments.

5.1 Polarisation

- Definition of OFF/ON Values for Power Supplies

DEFOFF {<module> = <off position>}

or

DEFON {<module> = <on position>}

DEFON HFI1=1.5 HFI2=1.0 HFI3=1.5

- Definition of Spin Turner IN and OUT (3D polarisation, only E1)

DEFTURN YIN = <9 values> *ZIN* = <9 values> *VFI1* = <9 values> -

YOUT = <9 values> *ZOUT* = <9 values> *HFI1* = <9 values>

- Definition of a Polarisation Cycle with POLB (polarisation begin) and POLE (polarisation end)

POLB

SET {<polarisation_module> = <target position>} if V6

SET {<polarisation_module> = ON/OFF} otherwise

MES [<presetvalue>]

...

POLE

POLB

SET HF=OFF VF=OFF FL1=OFF FL2=OFF

MES

SET HF= ON VF=OFF FL1=OFF FL2=OFF

MES 2500

POLE

- Switch ON/OFF Polarisation to Execute the Polarisation Cycle

POL ON/OFF

POL ON

Note:

The definition and execution of a measuring cycle (POLB, POLE, POL) may be freely used (independent of polarisation!).

5.2 Coordinate Transformations (Angles \Leftrightarrow H,K,L, ...)

- Lattice Parameter

LAT <6 values>

LAT 7.235 7.235 1.8242 90 90 90

- UB matrix

(number of significant parameters depends on instrument command e.g. E2 with FLAT: only 1st, 2nd, 4th, 5th parameter meaningful)

UBM <9 values>

UBM ...

- Wave Length

WAV <wavelength>

WAV 2.43138

- Method of Conversion H,K,L,E into Angle Values

(conversion method depends on instrument and EXPTYPE command)

MSH <conversion method>[, <alternative>]

with conversion method:

0: bisectric (PHIS \neq 0), or 4 circles

1: bisectric (PHIS = 0), or 2 circles

2: parallel

3: multi counter inclination (flat cone)

4: 2 circles, *OMGX*

5: goniometer (only E4)

MSH 0

MSH 1, 2¹

- Passive Modules (angle positions calculated from e.g. H, K, L will not be driven)

PASS {<module>}

PASS TTHS

- Calculate Angles

CALC {<module> = <value>}

¹alternative 2 is selected

CALC H=1 K=2 L=1 E=0.1

- Calculate H, K, L (only for A2, EULER)

CALCHKL OMGS = $\langle value \rangle$ TTHS = $\langle value \rangle$ PHIS = $\langle value \rangle$ CHIS = $\langle value \rangle$

CALCHKL OMGS=12.4 TTHS=2.2 PHIS=16.8 CHIS=0.1

5.3 Diverse special

- D Spacing monochromator[, analyzer]

DS $\langle value \rangle$ [$\langle value \rangle$] different number of significant parameters

DS 3.3542 3.3542

- Scattering Sense monochromator, sample[, analyzer]

1: left, -1: right (reverse for E10)

SS $\langle 2 values \rangle$ [$\langle value \rangle$] different number of significant parameters!

SS 1 1 1

- Fix initial (0) or final (1) Energy wave vector, wave length

FX 0/1

FX 1

- Energy Unit

(0: THz, 1: meV)

EU 0/1

EU 1

- Definition of Resolution Steps

RESOSTEP $\langle number of steps \rangle$

RESOSTEP 15

- Definition of Partial Resolution Acquisition

RESOPART $\langle reso_steps \rangle \langle start_step \rangle \langle end_step \rangle$
 $(start_step \langle = end_step \langle = reso_steps)$

RESOPART 25 13 25

- Select the omega module used for scans and conversions

OMEGA_MOD $\langle omega module \rangle$

OMEGA_MOD OMEGAX

Note:

The default is *OMGS*; a definition could be reset via *RESET OMEGA_MOD*.

- Coupling set modules

CSM {*modules*}

CSM *TTHS* *FKM*

CSM *TEMP* *MAGF* *POSX*

Note:

This command defines dependencies of modules according to a lookup table; multiple sets of different modules are allowed; subsequent set and scan commands will use these definitions to determine the setpoint positions of dependend modules; the set command must not contain more than one module of a csm definition; all sets are cleared via *RESET CSM*.

- Select action after a fatal error (a device requested a RESET)

FATAL_ERROR *<action>*

with action:

IGNORE: try to continue scan

HALT: (default) halt scan and ask user

ABORT: abort scan and wait for next command

FATAL_ERROR *HALT*

6 Data acquisition and other actions

This set of commands causes actions with real hardware access, e.g. readout or movement to new set points of devices.

6.1 Perform Basic Actions

- Read and Print Current Positions

READ {<module>}

READ OMGS TTHM TEMP

- Drive Modules to Target Position

SET {<module> = <value or identifier>}

SETC {<module> = <value or identifier>} *continue / no waiting*

SET OMGS=3.23 PHIS=-23.02

SET H=1 K=2 L=1

SET HF=ON, MATRIX=XY

- Data Acquisition

MES [<preset>]

MES

MES 2500

Note:

This command starts a single data acquisition without incrementing the run number and overwrites the temporary measurement file! For regular measurement files refer to (next) section scans.

A given preset is temporary used instead of the default defined by command MM1.

6.2 Scans

- Scan around a central position

SGEN STEP = <steps>{<module> = <central_pos>}{<d_module> = <step_width>}

SGEN STEP = 11 OMGS=10.3 DOMGS=0.2 TTHS=20.6 DTTHS=0.4

```
SGEN STEP = 21 DOMGS=0.21
SGEN STEP=31 H=1 K=2 L=1 DH=0.1 DK=0 DL=0
```

- Scan from start to end position

```
SGEN1 STEP = <steps>{<module> = <start_pos><end_pos>}
```

```
SGEN1 STEP =21 OMGS=20.5, 31.8
SGEN1 STEP=13 H=1, 1.2 K=2,2.2 L=3,3.2 E=0
```

- Definition of an arbitrary Scan using SCANB (scan begin) and SCANE (scan end)

```
SCANB {<action>} SCANE
```

```
SCANB
SET OMGS=5.3
SET MAGF=2.5
MES
```

```
...
SCANE
```

Note:

It is not allowed to change (redefine) the master module (MM1) inside the scan.

Continuous scans There are two variants of continuous scans. Both at first set *<motor>* to the position *<start>* and subsequently continue to move *<motor>* to position *<target>*. At the moment *<motor>* reaches position *<start_mes>* CARESS starts the measurement. The scan terminates if *<motor>* reaches position *<target>* or stops for any other reason.

- Continuous Scan with *fixed* number of steps

```
SGEN2 <motor> = <start>, <start_mes>, <end_mes>, <target> STEP = <steps>
```

```
SGEN2 OMGX=11.7, 12.1, 42.8, 44.3 STEP=27
```

Note:

There is no master counter involved (MM1 is to be ignored!) and all counters are used as slave modules. CARESS is doing *<steps>* measurement steps without stopping the movement nor the counters. The parameter *STEP* is defined to be the last parameter of the command.

- Continuous Scan with *variable* number of steps

```
MM1 <master counter> = <preset>
SGEN2 <motor> = <start>, <start_mes>, <end_mes>, <target>
```

```
MM1 TIM1=10000
SGEN2 OMGX=11.7, 12.1, 42.8, 44.3
```

¹scan around actual position of OMGS

Note:

The master counter (MM1) is used to determine the measurement steps while all slave counters and the *<motor>* are continuously running.

6.3 Initialisation

- Initialise Hardware

INITH

INITH

Note:

This command updates all CARESS internal information about the hardware configuration. Use *ini th* to read in updated hardware description files.

- Reset / Delete Command Parameters

with commands: CALIB, CSM, OFL, ONL, PASS, PERMOVE, PRO, PRW, PRWA, RESOSTEP, SCANV, SELDET, SL1, SPEC, TOF, WAIT1, WAIT2, WIND, WINDA

RESET {*<commandkey>*}

RESET PASS

7 Special actions

This section lists commands to deploy some special hardware setups/access.

7.1 Reference movements

- Move one or more motor axis, each to one reference point

REF {*module*} = *direction*}
with *direction=1* \equiv *positive limit switch*,
direction=2 \equiv *negative limit switch*

REF OMEGA=1, SLIT1_X=2

- Move one or more motor axis, each to one reference point and define this point as given position value

REFSET {*module*} = *direction*, *position*}
with *direction=1* \equiv *positive limit switch*,
direction=2 \equiv *negative limit switch*,
position=new position value after movement

REFSET OMGS=1, 27
REFSET OMGS=2, 127 TTHS=1, -21.7

- Definition of a permanent 'Movement' of a magnetic field

PMOVE *module* = *ccu**ccd**wcu**wcd**rcu**rcd*
with: *ccu/ccd* = *cycle current up/down*
wcu/wcd = *wait current up/down*
rcu/rcd = *rate current up/down*

PMOVE MAGF=0.5 1 0.7 1.22 2

Note:

This command currently works for magnetic fields controlled by a bruker bec-1 power supply (CARESS device kind 68).

7.2 Universal device access

- receive an ASCII string from a gpib (iec bus) device

READGPIB $\langle X \rangle : \langle Y \rangle$

with X = logical unit number of the gpib controller and
 Y = address of the gpib device

or

READGPIB $\langle Z \rangle$

with Z = logical unit number of the gpib device.

READGPIB 1:30

READGPIB 27

- send an ASCII string to a gpib (iec bus) device

WRITEGPIB $\langle X \rangle : \langle Y \rangle$ 'string'

with X = logical unit number of the gpib controller and
 Y = address of the gpib device

or

WRITEGPIB $\langle Z \rangle$ 'string'

with Z = logical unit number of the gpib device.

WRITEGPIB 1:19 "MEASure:VOLTage:DC? 20, 0.001"

WRITEGPIB 27 "SETup:TXPower:COUNT:STATE ON"

- send an ASCII string to a device/module

LOADTEXT $\langle module \rangle = \langle string \rangle$

$\langle string \rangle ::= \langle identifier \rangle | \langle text \rangle$ (no ' in text)

LOADTEXT ADET='LISTMODE ON'

Note:

To interpret the sent information is exclusively up to the device (driver).

- send the contents of a file to a device/module

LOADFILE $\langle module \rangle = \langle filename \rangle$

$\langle string \rangle ::= \langle identifier \rangle | \langle filename \rangle$

LOADFILE ADET='adet-special.cfg'

- Set an input register bit to either logical 0 or 1

SETBIT $\langle bit - number \rangle \langle value \rangle$

SETBIT 2 1

Note:

This command is used at a special wired register at V4. (cf. WAIT4BIT at page 24)

7.3 Complex actions

- Sophisticated command at instrument M2

SETCM2 <device><motoraxis><setpoint><rangelow><rangehigh><mode><interval>

SETCM2 COU, TTHS, 123, 5, 124, 0, 1.7

Note:

This command initiates a very special algorithm at instrument M2 by sending several data to the named device. It starts a measurement within range while moving from current to set point position of motor axis. mode and interval are additional measurement parameters.

This command is blocking (see ENDCM2 at on this page)

- Terminate the SETCM2 command

ENDCM2

ENDCM2

Note:

Since SETCM2 is a blocking command it is necessary to terminate it explicitly with ENDCM2.

8 Control special hardware devices

A set of commands to deal with more sophisticated hardware systems is shown in the following paragraphs.

8.1 Definition of Compound Positions, etc.

- Sample Changer Positions (V4)

```
DEFPOS {<pos_no><pos_x><pos_y>}
```

```
DEFPOS 1 2.3 0, 2 0.5 0, 3 0.52 0.1, ...
```

- Define List Position

```
DEFLPOS LPOS = <pos_no>{<module> = <value>}
```

```
DEFLPOS LPOS=1 H1=2.0 H3=1.5 H10=-0.5 ...
```

- Change (defined) List Position

```
CHGLPOS LPOS = <pos_no>{<module> = <value>}
```

```
CHGLPOS LPOS=1 H3=1.6
```

- Define and name special Positions (E2, E7, V4, M1, M2) or Distances (V6)

```
DEFINE {<module>_<name><value>}
```

```
DEFINE BSX_OUTX 53.1 BSY_OUTY 30.95
```

8.2 Define Data Acquisition Unit DAU

The DAU is a special hardware developed at Hahn-Meitner-Institut Berlin and used at instruments V5/SPAN only.

- Enable detectors and define data width (16 bit or 32 bit) for incrementing

SI: single detector banks and/or

MU: area detector,

I2: 16 bit (2 bytes), I4: 32 bit (4 bytes)

```
CONFIG {<detector> = <width>}
```

```
CONFIG SI=I4 MU=I2
```

```
CONFIG SI=I2
```

- Define Time-Of-Flight Parameters (note the DAU has 384KB memory for detector data)

TB: time base (0..4: 0=125ns/8MHz, 1=250ns, 2=0.5 μ s, 3=1 μ s, 4=2 μ s),

DT: delay time (units of TB for delay between pick-up signal and start of time window),

CN: number of channels (1..4096),

CW: channel width in units of TB ($\langle \rangle$ =1)

TOF TB = $\langle value \rangle$ DT = $\langle value \rangle$ CN = $\langle value \rangle$ CW = $\langle value \rangle$

TOF TB=0 DT=20 CN=1024 CW=10

defines TOF mode and TOF parameters

time base 125ns (8MHz), delay time 2.5ms, 1024 channels of 1.25ms

- Switch from TOF Mode to NSE Mode

RESET TOF

RESET TOF

8.3 Neutron Chopper System

- This command is used to set up several parameters of a chopper system. The real activation of these parameters is done with 'SET CHOP = RUN'

Define State

0: synch,

1: asynch,

Speed (RPM),

Speed Ratio and

Phase in relation to Chopper 1

CHOPPER STATE = {{ $\langle chop_no \rangle$ $\langle state \rangle$ } SPEED = {{ $\langle chop_no \rangle$ $\langle speed \rangle$ }

RATIO = {{ $\langle chop_no \rangle$ $\langle ratio \rangle$ } PHASE = {{ $\langle chop_no \rangle$ $\langle phase \rangle$ }

CHOPPER STATE=2 0, 3 0, 4 0, SPEED=1 8000 -

RATIO=2 4 PHASE= 2 0, 3 0, 4 0

9 Flow control

A group of commands to control the actual command execution is to find here.

9.1 Wait commands

- Wait for seconds (unconditional)

WAIT $\langle seconds \rangle$

WAIT 10

- Wait before start of a scan

WAIT1 $\langle seconds \rangle$

WAIT1 20

- Wait between accomplished drive and start of acquisition in a scan (repeatedly)

WAIT2 $\langle seconds \rangle$

WAIT2 15

- Wait for an external signal (hardware, i/o-module trigger)

WAITTRIG

WAITTRIG

- Wait for either logical 0 or 1 at an input register bit

WAIT4BIT $\langle bit - number \rangle \langle value \rangle$

WAIT4BIT 2 1

Note:

This command is used at a special wired register at V4. (cf. SETBIT at page 20)

- Wait until logical expression becomes true

WAIT4EXPR $\langle expression \rangle$

WAIT4EXPR ABS(TEMP-T_SAM) \leq 0.25

9.2 Input Utilities

Used for Variables, Procedures, Loops, Dummy Command Parameters

- Definition of Variables, Arithmetic Operations (+, -, *, /), Mathematical Functions (sin, cos, tan, asin, acos, atan, sinh, cosh, tanh, exp, log, log10, ceil, floor, sqrt, fabs) and Assignment

LET $\langle \text{variable} \rangle = \langle \text{expression} \rangle$

```
LET A = 5.4
LET START = -1
LET WIDTH = (2+OMGS)*SQRT(A) + START
```

- Definition of Procedure

PROCEDURE $\langle \text{proc_name} \rangle (\langle \text{formal_parameter_list} \rangle)$
sequence of CARESS commands
ENDP

```
PROCEDURE MYPROC(STEPS, CENTRE, STEP_WIDTH)
SGEN STEP=STEPS OMGS=CENTRE DOMGS=STEP_WIDTH
...
ENDP
```

- Call of Procedure

DO $\langle \text{proc_name} \rangle (\langle \text{actual_parameter_list} \rangle)$

```
DO MYPROC(11, 5, 0.1)
```

- List currently defined Assignments or Procedures

LIST $\langle \text{command} \rangle$ (LET or PROCEDURE)

```
LIST LET
```

- Loops (loops may be nested)

REPEAT $\langle \text{repeat count} \rangle$
{any sequence of commands}
ENDR

```
LET I = 0
LET H1_START = 1.5
LET H1_DELTA = 0.2
REPEAT 10
LET H1_POS = H1_START + I*H1_DELTA
SET H1 = H1_POS
MES
LET I = I+1
ENDR
```

- Open a Command Parameter File
PARF <filename> !(max.20 characters, CARESS appends extension .par)
 PARF myhkl
- Replace dummy Command Parameters (\$) by Values from the Command Parameter File myhkl.par
SET commandparameters
 SET H=\$ K=\$ L=\$
 SET FL=\$
- Exit active Loop, if end-of-record (EOR) or end-of-file EOF
IF <flag> EXIT
 REPEAT 10
 SET H=\$ K=\$ L=\$
 MES
 IF EOR EXIT
 IF EOF EXIT
 ENDR

9.3 Control CARESS by an external user program, ...

- Call External User Program
EXTERN <text>
 EXTERN myprog
 EXTERN center
- Call External PV-WAVE User Program
EXTERN WAVE <text>
 start PV-WAVE program myprog.pro:
 EXTERN WAVE myprog
 start PV-WAVE (and type in PV-Wave commands):
 EXTERN WAVE
- Return from external control to normal CARESS control
INTERN
 INTERN

Note:

INTERN is not allowed in command files, procedures, and repeat sequences.
 EXTERN is ignored, if CHECK is active.

EXTERN will interrupt an active command file.

If EXTERN becomes active in a procedure or repeat sequence out of a command file, no command file may be opened by the external process.

An external process may not define/use procedures nor repeat sequences (it will be more convenient to use the facilities of the programming language!).

- Set a semaphore and abandon current control

BLOCK <GID> (with $0 \leq GID \leq 9$)

BLOCK 1

Note:

This command is used for synchronisation reasons in a multi-client setup. BLOCK abandons the current control and defines a semaphore group. The **only** way to get control back is to use UNBLOCK (description on the current page) from another CARESS client.

- Clear semaphore in multi-client environments

UNBLOCK [GID]

UNBLOCK

UNBLOCK 3

Note:

This command is used for synchronisation reasons in a multi-client setup. UNBLOCK clears all semaphores, with an optional group id only this semaphore group is cleared. Semaphores are defined with BLOCK (description on this page).

9.4 Interrupt / Kill / Continue / Terminate CARESS Activities

- kill current activit

KILS

KILS

- delete data file (after KILS, only)

KILD

KILD

- cancel command input queue including all open command file

KILC

KILC

- cancel active repeat loops and/or procedure

KILP

KILP

- continue normal processing

CON

CON

- abort a cancel action (KILS, KILC, KILP, INTERN and CON) (command version, only)

ABO *<killcommand>*

ABO KILP

- abort a cancel action (KILS, KILD, KILC, KILP, INTERN and CON)

ABOD *<killcommand>*

ABOD KILD

- terminate all CARESS activities (not allowed in command files, procedures, and repeat sequences)

EX or *EXIT*

EXIT

Note:

To interrupt an activity click the HALT button within the GUI or type CNTRL/h within the command version of CARESS.

10 Miscellaneous commands

Some Commands do not fit into categories before and will be listed here.

- Check commands

CHECK ON/OFF [options]

CHECK ON
@mycommands
CHECK OFF

Options to restrict check depth are:

MAXREPEAT = <n> limit repetition of loops to <n>

range (aliases to *MAXREPEATS* and *MAXLOOP*)

MAXFILE = <n>

range limit nesting script level to <n>

range (aliases to *MAXFILES* and *NESTED*)

NEWONLY = ON/OFF

range if *ON*, include script files only once

range (aliases to *NEWFILEONLY*, *NEWFILES ONLY*

range and *UNIQUE*)

Note:

The enclosed commands are only checked, the parameters are not stored! Procedures are being defined, independent of CHECK! The enclosed commands and CHECK OFF are not written into the CARESS session protocol. CHECK ON/OFF is not allowed in command files, procedures, and repeat sequences. Typing CHECK OFF within procedure or repeat definition and CHECK active will cancel all open procedure/repeats. A warning message is printed. By this mechanism errors in command files (e.g. missing ENDR) may be detected, if checked by the GUI. Caution! Don't forget to edit the faulty command file before executing it!

- Switch OFF/ON Archiving of Data on central archive system

ACV OFF/ON

ACV OFF

- Bypass Condition(s) which prevents Data Acquisition e.g. SHUTTER closed (V4)

BYPASS ON|OFF

BYPASS ON

- Disable *keep position* of a Motor Control or all Motors
i.e. allow persistent manual positioning (Y66 motor control)

MANUAL [⟨module⟩]

MANUAL

MANUAL OMGS

- Set Modules "offline"

OFL {⟨module⟩}

OFL TTHA

- Disable/enable Protected Input
(CARESS command version only, not for GUI)

PASSW ON/OFF

PASSW ON

- List the current Occurrence of a Command

SHOW ⟨command key⟩

SHOW BLS

Index

ABO, 28
ABOD, 28
ACQPMOVE, 8
ACV, 29

BLOCK, 27
BLS, 7
BYPASS, 30

CALC, 13
CALCHKL, 14
CALIB, 8
CHECK, 29
CHGLPOS, 22
CHOPPER, 23
COM, 7
CON, 28
CONFIG , 22
CSM, 15

DEFINE, 22
DEFLPOS, 22
DEFOFF, 12
DEFON, 12
DEFPOS, 22
DEFSAM, 6
DEFTURN, 12
DISP, 10
DO, 25
DS, 14

ENDCM2, 21
ENDP, 25
ENDR, 25
EU, 14
EX, 28
EXIT, 26, 28

EXPTYPE, 6
EXTERN, 26

FATAL_ERROR, 15
FX, 14

IF, 26
INITH, 18
INTERN, 26

KILC, 27
KILD, 27
KILP, 28
KILS, 27

LAT, 13
LET, 25
LIST, 25
LOADFILE, 20
LOADTEXT, 20
LOGDEV, 10

MANUAL, 30
MES, 16
MM, 7
MM1, 7
MRELA, 9
MSH, 13

OFL, 30
OMEGA_MOD, 14
ONL, 6

PARF, 26
PASS, 13
PASSW, 30
PEAK, 10
PMOVE, 19

POL, 12
POLB, 12
POLE, 12
PRO, 7
PROCEDURE, 25
PRT, 9
PRTMOD, 9
PRW, 10
PRWA, 10

READ, 16
READGPIB, 19
REF, 19
REFSET, 19
RELA, 9
REPEAT, 25
RESET, 10, 18, 23
RESOPART, 14
RESOSTEP, 14

SAM, 6
SCANB, 17
SCANE, 17
SCANV, 8
SELDET, 8
SET, 16
SETBIT, 20
SETCM2, 21
SGEN, 16
SGEN1, 17
SGEN2 (fix.), 17
SGEN2 (var.), 17
SHOW, 30
SL1, 7
SOF, 7
SPEC, 10
SS, 14

TITLE, 6
TMO, 7
TOF, 23
TOL, 7

UBM, 13
UNBLOCK, 27

USN, 6

WAIT, 24
WAIT1, 24
WAIT2, 24
WAIT4BIT, 24
WAIT4EXPR, 24
WAITTRIG, 24
WAV, 13
WIND, 8
WINDA, 8
WRITEGPIB, 20

User's Notes