



# SPEC - REFLECTOMETER CONTROL



# SPEC - REFLECTOMETER CONTROL

Last change: 2011/05/17 AG/FS

Contacts: Franz Schäfers (12946), Ingo Packe (12943), Andreas Gaupp (12957)

## 1. Getting started

- Log-in on POLARIS.EXP.BESSY / UNIX-system as: *SPECUSR*
- Password: *brillanz*
- Open a terminal-window
- Type: *spec*
- Enter two character for your user name (case sensitive)
- To modify volume and/or file number use “*Volume*” and “*Filenumber*”
- **SPEC and UNIX are case-sensitive!!!**

## 2. Available Devices

<b>Motors:</b>	<b>function/device</b>	<b>mnemonic name</b>	<b>(unit)</b>	<b>hardware</b>	<b>soft limits</b>
1.	theta	tha	(deg.)	EPICS REFLC:h0000007 / VME #8	-20 - 140
2.	twotheta	twt	(deg.)	EPICS REFLC:h0000006 / VME #7	-45 - 230
3.	detector	det	(deg.)	EPICS REFLC:h0000005 / VME #6	-45 - 15
4.	sample	smp	(mm)	EPICS REFLC:h0000004 / VME #5	1 - 49
5.	filter	flt	(mm)	EPICS REFLC:h0000003 / VME #4	1 - 140
6.	mirror	mirror	( $\mu$ rad)	EPICS Macro-Mot / PM4M1:RxAbs	-50 - 50
7.	mono	mono	(eV)	EPICS / Macro-Mot d0812pgm1:mono...	10 - 5e6
8.	V1	V1	(Volt)	Keithley #10 Voltage source Macro-Mot	-100- 100

### Counters:

1.	Io Au-mesh	kth1	(A)	GPIO / Keithley 617 #10
2.	GaAsP-diode	kth2	(A)	GPIO / Keithley 617 #11
3.	ring current	rc	(mA)	EPICS
4.	Channeltron	counter1	(cts)	counter card
5.	Counter	counter2	(cts)	counter card

## 3. Device handling: move motors, read detectors

If nothing moves, type: *reconfig*

**mtha, mtwt, mdet, msmp, mflt, mmono, mmirror, mV1**

Description: Moves motor to destination  
Usage: e.g. *mmono destination (eV)*

**mrtha, mrtwt, mrdet, mrsmp, mrflt, mrmono, mrmirror**

Description: Moves motor relative to current position  
Usage: e.g. *mrmono eV*

**wa**  
Description: Shows the actual positions of all motors in user and dial coordinates  
Usage: *wa*

**read\_keithley**  
Description: Returns the actual value read from keithley  
Usage: *read\_keithley kthl*

**ct**  
Description: Reads all counter channels  
Usage: *ct*

## 4. Scan macros

**bragg**  
Description: Performs a Bragg scan of  $\theta$  in reflection. Scans around the Bragg angle by  $\pm\Delta\theta$  for energies in given range; d-spacing in nm.  
Usage: *bragg delta\_theta steps\_theta start\_energy stop\_energy steps d-spacing*

**detector**  
Description: Performs a detector scan  
Usage: *detector start stop steps*

**energy**  
Description: Performs an energy scan  
Usage: *energy start stop steps*

**Energy**  
Description: Performs an energy scan. Prior to scan user must determine proper position of mirror at both ends of scan  
Usage: *energy start stop mirror\_start mirror\_stop steps*

**filter**  
Description: Performs a scan of filter-feedthrough to find the filter positions  
Usage: *filter start stop steps*

**gr\_onblaze**  
Description: Performs an on-blaze energy scan of a grating at a certain order  
Usage: *gr\_onblaze energy\_start energy\_stop steps l/mm blazeangle order*

**mirrorRx**  
Description: Performs a scan of Rx (tilt, incidence angle) of premirror M1  
Usage: *mirrorRx start stop steps*

**ml\_peakref**  
Description: Performs an energy scan on the Bragg maximum of multilayer  
Usage: *ml\_peakref energy\_start energy\_stop steps d-spacing (nm)*

**sample\_x**  
Description: Performs a sample scan in x (perpendicular to plane of reflection)

Usage: *sample\_x start stop steps*

### **stripchart**

Description: Mimics a stripchart recorder, I.e. performs a time scan  
Usage: *stripchart timeperstep steps*

### **theta**

Description: Performs a theta scan  
Usage: *theta start stop steps*

### **twotheta**

Description: Performs a two-theta scan  
Usage: *twotheta start stop steps*

### **twotoone**

Description: Performs a 2:1 scan that scans tha and twt with twt moving twice the range of tha  
Usage: *twotoone theta\_start stop steps*

### **TwoToOne**

Description: Performs a 2:1 scan that scans tha, twt and det with twt moving twice the range of tha. Prior to scan user must determine proper position of detector at both ends of scan  
Usage: *TwoToOne theta\_start stop steps det\_start stop*

### **Vscan1**

Description: Performs a voltage scan with Keithley #10  
Usage: *Vscan1 V\_start V\_stop steps*

## **5. Miscellaneous – SPEC-symbols, commands, macros (a selection)**

### **abort, stop, resume a scan**

Usage: *p* then “*r*” (resume) or “*ctrl/c*”

### **autosave**

switches autosave feature on or off  
Usage: *autosave=“on”* or *autosave=“off”*

### **check\_beamshutter**

starts scan only when beamshutter is open (set by default)

### **Comment**

Shows and sets the comment for the datafile

### **config**

only for experts, must be SPECADM  
Usage: *config*  
to quit: *Cntr C*, to write: *w*, to edit: “*”*, to toggle: *c*

### **dweltime**

Shows and sets the counter dwelltime, Keithley preamp settling time and sequence of Keithley readout.  
We wait PreampSettlingTime before starting the loop for number of Readouts. The loop over Readouts also waits for PreampSettlingTime before taking any reading.

		Keithley specification: 2.5 sec in 20 pA and 200 pA 15 msec in 2 nA, 20 nA, 200 nA 2 msec in microA and mA ranges
	Usage:	<i>dwelltime Dwelltime PreampSettlingTime</i>
<b>delay</b>		Shows and sets the mechanical delaytime (sleep) in seconds before counters are read
	Usage:	<i>delay</i> or <i>Delay delaytime</i>
<b>Detector</b>		gives info on available detectors (channel 2) and their relative position
<b>filenumber</b>		Shows and sets file number for current user
	Usage:	<i>filenumber</i> or <i>Filenumber</i>
<b>Filterinuse</b>		moves one out of 6 filters into the light beam
<b>help</b>		Invokes the SPEC on-line help
<b>home</b>		
	Description:	clears the lightpass in reflectometer for the SURICAT experiment (tha: 0, twt: 0, det: -15, smp: 49)
	Usage:	<i>home</i> or <i>reflec_home</i>
<b>ignore_beamshutter</b>		starts scan without checking beamshutter
<b>Io</b>		
	Description:	moves sample out of and detector into the light beam for Io measurement (smp: 49, tha: 0, twt: 0, det: 0)
	Usage:	<i>Io</i>
<b>pinhole</b>		gives info on available pinholes and their position on manual feedthrough
<b>plot</b>		Plots the last scan
<b>plotsselect</b>		Selects the channels to plot on line (max 3)
	Usage:	<i>plotsselect counters</i>
<b>prdef</b>		prints definition of macro
<b>qdo</b>		Includes macros into SPEC
	Usage:	<i>qdo usermacros.mac</i> or <i>qdo spec.mac</i>
<b>quit</b>		EXIT from the program
<b>Readouts</b>		Shows and sets the number of readouts, (de)activates Keithley, and determines readout sequence for ReadoutSequence = "sim" : reading pseudo simultaneously, ReadoutSequence = "seq" : reading sequentially (only effective for Readouts > 1). "sim" recommended since PreampSettlingTime can be reduced by time it takes to read the Keithleys. "seq" recommended when digitising noise etc. is dominating.

	Usage:	<i>readouts</i> or <i>Readouts</i> <i>ReadoutSequence</i> = <i>sequential</i> or <i>simultaneous</i>
<b>Sample or sample</b>		Shows and sets the sample comment/description for the datafile
<b>Sample_info</b>		Shows the motor-positions on sample holder
<b>savedata</b>		Saves the data of the last scan to POLARIS\DATA-directory (no overwriting of existing data files), auto increment of Filenumber.
	Filename:	username_ref_volume.filenumber (e.g.: fs_ref_0.1)
	Usage:	<i>savedata</i> or <i>save_data</i> or <i>datasave</i> or <i>data_save</i>
<b>set</b>		Redefines a motorposition in user coordinates
	Usage:	<i>set motor position</i>
<b>set_dial</b>		experts only
	Usage:	<i>set_dial motor position</i>
<b>setplot</b>		Changes plot options. Attention: Default-setting is NO PLOT and NO HIGH RESOLUTION PLOT!!!
	Usage:	<i>setplot</i>
<b>syms</b>		displays all SPEC symbols, arrays, strings,...
<b>TAB</b>		displays all 641 SPEC commands
<b>volume</b>		Shows and sets volume number for current user
	Usage:	<i>volume</i> or <i>Volume</i>
<b>whoami</b>		Shows present user name
<b>File usermacros.mac</b>		write your own macro-routine for your measurement sequence into "usermacro.mac" e.g. according to: <pre>def refler_home '{ # moves reflectometer sample out of light beam # to give beam to SURICAT exp. station mtha 0 mtwt 0 mdet -15 msmp 49 }'</pre>
	Include macro:	<i>qdo usermacros.mac</i> or <i>qdo spec.mac</i>

## 6. Trouble shooting

**reconfig**            updates current configuration (similar to new start)

## 7. Data transfer --- Connection to BESSY network drives

### Available Networks in BESSY Experimental hall

*USR* (yellow label) for User computers (100 Mbit)  
*EXP* (white label) for BESSY devices (100 Mbit)  
*BLC* (generally not accessible)  
*WLAN* (10 Mbit, at request only)

### Login Username / Password

*monop / monop*  
*DIP\_KMC1\$ / DIP\_KMC1*  
or: *bessyguest / bessyguest*  
or: *teas\yourname / yourpassword*

### Available data-drives

<a href="\\fs.exp.bessy.de\scratch">\\fs.exp.bessy.de\scratch</a>	(user/password: <i>monop/monop</i> )
<a href="\\fs.exp.bessy.de\public">\\fs.exp.bessy.de\public</a>	( <i>monop/monop</i> )
<a href="\\bessy.exp.bessy.de\monop">\\bessy.exp.bessy.de\monop</a>	( <i>monop/monop</i> )
<a href="\\fs.exp.bessy.de\linse">\\fs.exp.bessy.de\linse</a>	( <i>yourname/passowrd</i> )
<a href="\\polaris.exp.bessy.de\specusr">\\polaris.exp.bessy.de\specusr</a>	( <i>specusr/brillanz</i> (Reflectom.))
<a href="\\speckle.exp.bessy.de\specusr">\\speckle.exp.bessy.de\specusr</a>	( <i>specusr/brillanz</i> (Polarimeter))

etc...  
no connection to OS2-drives possible

### Connection to network from

<b>Windows</b>	Open	<i>Arbeitsplatz – Extras – Netzlaufwerk verbinden</i> <i>Your workspace – Extras – Map network drive</i>
<b>Windows or OS2</b>	Open	<i>Eingabeaufforderung</i>
	Check online status:	<i>ping bessy (fs, polaris, speckle...)</i>
	Check network drives:	<i>net use</i>
	Check IP-configuration:	<i>ipconfig</i>
	Create new drive (e.g.):	<i>net use X: <a href="\\bessy\monop">\\bessy\monop</a></i>
<b>UNIX</b>	Open	<i>HOME window</i>
	Enter location e.g.:	<i>smb://bessy/monop</i>
	Enter username e.g.:	<i>teas\monop</i>
	Password e.g.:	<i>monop</i>

## 8. Connection to BESSY printers

**Printer server:** <\\diprint\printername> (user: *teas\monop / monop*)  
**Windows:** *START – run – <\\diprint> – user/password – select printer from list*

e.g. printer in the upper level next to KMC-1: **waz\_ps** or **waz\_pcl6**  
UE56/2 PGM1,2: **nnn\_ps** or **nnn\_pcl6**  
Reflectometer: **nkr.prn.bessy.de**

## 9. PM-4 Optics Beamline settings

### for SURICAT:

<b>Apertures:</b>	4.8/2.0/7.0/16.0
<b>Mirror M1 Rx -Motor 5:</b>	0 +/- 15 $\mu$ rad
<b>Monochromator:</b>	360 l/mm ((Pos. 74 mm on hand-operated LMD)
$c_{ff}$ :	3
<b>Reflectometer:</b>	SPEC-control: reflc_home (theta=0°, twotheta=0°, detector=-15°, sample_x=49 mm, filter=10 mm)
Pinhole feedthrough:	50 mm
Io-mesh feedthrough:	0 mm
Filter feedthrough:	10 mm
Vacuum:	$<2 \times 10^{-6}$ mbar

### for REFLECTOMETER

<b>Front end Apertures:</b>	2.3/-1.0/6.5/12.5 mm (=3.3 x 19 mm <sup>2</sup> (v x h) @14.6 m = 0.23 x 1.3 mrad <sup>2</sup> ) (Offset-values:-4.5/-4.9/-3.9/-5.6)
<b>Mirror M1 Rx - Motor 5:</b>	0 $\mu$ rad +/-15 $\mu$ rad (energy-dependent)
<b>Straylight apertures</b> (25 mm hand-operated LMDs)	14/17 mm (top/bottom) 8 / 4 mm (left/right)
<b>Monochromator:</b>	1228 l/mm (Pos. 22 mm on hand-operated LMD) $c_{ff}$ : (recommended) $\leq 2$
<b>Straylight apertures</b> (25 mm hand-operated LMDs)	15/15 mm (top/bottom)
<b>Intermediate focus apertures</b> (25 mm hand-operated LMDs)	14/12 mm (top/bottom) 18/12 mm (left/right)

## 10. Sample change – Venting of Reflectometer

1. Move all motors to *HOME* position by SPEC-program. **Especially the FILTER-motor must be at home-position 10 mm for protection.**

2. Close hand-valves on both sides of reflectometer (to beamline and to SURICAT).

**Note: both valves are double valves, close the inner ones, next to the chamber (labeled with “Fenster MgF2” and “Fenster”, resp.) !!!**

3. Press the “STOP”-button on the turbopump control (Turbotronik NT 340M).

**Note:** the valve (NW150 CF) between turbopump and chamber closes automatically.

4. Stop the roughing pump unit (DCU Pfeiffer Vacuum) on the other side of the chamber.

5. Switch off the Penning pressure gauge B1 (at BALZERS TPG 300 press buttons “Sensor” and “Step”). Change pressure reading to the pirani gauge at B2 (press button “Sensor”).

**Now you are ready to vent the chamber via the nozzle-valve with gaseous Nitrogen from the big bottle next to the chamber.**

6. Open the valve on top of the N<sub>2</sub>-bottle.

7. Open the nozzle-valve on top of the reflectometer chamber slowly while checking the pressure in the chamber. It should increase up to approximately 700 mbar (is not calibrated).

**Note:** check vacuum of last beamline section in the yellow beamline rack (Ion pump power supply labelled “Austrittspalt”). Vacuum should stay in 10<sup>-9</sup> mbar range. If not, you have a leak to the beamline. stop venting immediately, close the nozzle valve.

8. Release latches at the load-lock door (CF 150 Window) on the chamber. At atmospheric pressure the door opens automatically.

9. Close the valve on top of the N<sub>2</sub>-bottle. Close the nozzle-valve at the chamber.

**Now you can change sample**

## 11. Sample change – Pump down of Reflectometer

1. Close nozzle valve tightly. Close load-lock door and fix the latches tightly.
2. Start both the roughing pump unit (DCU Pfeiffer) and the Scroll-roughing pump (VARIAN). Change the 3-way cross between turbopump and roughing pump section to position 2 to pump down with the Scroll-pump.
3. Press the “START”-button on the turbopump control (TURBOTRONIK), then open the valve between turbopump and chamber by pressing the black knob on the valve control (Ventilsteuerung) next to the turbopump control.
4. Wait until the pressure reading B2 at BALZERS TPG 300 shows  $10^{-3}$  mbar and until the turbopump has come to full speed. If so, only three green lights remain shining on the control unit. This may take up to 15-30 minutes.
5. Check whether the roughing pump unit has come to full speed (1500 Hz). If so, change 3-way cross to position 1 to further pump with roughing pump unit. Switch off Scroll-pump.
6. Switch on pressure gauge B1 (on BALZERS TPG 300 press buttons “Sensor” and “func”). Should show  $10^{-4}$  mbar or better.
7. Wait 3 to 4 hours until vacuum is better than 5 times  $10^{-6}$  mbar. Now you can open the valve to the beamline.  
**Note: Check beamline vacuum (Austrittspalt). Should stay in the  $10^{-9}$  or  $10^{-8}$  mbar range.**

Now you can start to measure

