


DIP-09-2**KMC - 2**

CONTACT PERSON	building	room	phone	fax	e-mail
A. Erko	14.51	3371	2945	2990	erko@bessy.de
I. Packe	14.51	3354	2943	2990	packe@bessy.de

OPTICAL LAY-OUT (schematic)

KMC-2




KMC2.ds4 1.06.2001

premonochromator optics	no optics
entrance slit	Vertical and horizontal slits: 0-10 mm, water cooled
monochromator	Double-crystal monochromator, first crystal water cooled. 2 SiGe 111 graded crystals, angular range $2.16^\circ - 70^\circ$ Energy range 4 keV - 15 keV, feedback stabilization for EXAFS measurements
exit slit	Vertical and horizontal slits: 0-10 mm
postmonochromator optics	Rh-coated cylinder mirror with bending mechanics. Fixed horizontal focus at 35 m from the source, variable vertical focus 31 m - 36 m from the source. Vertical beam deviation: 0.5° above the horizontal plane
references	A. Erko, I. Packe, C. Hellwig, M. Fieber-Erdmann, O. Pawlitzki, M. Veldkamp, W. Gudat, "KMC-2: the new X-ray beamline at BESSY II, AIP Conference Proc. 521 , p.p. 415-418, (2000). A. Erko, I. Packe, W. Gudat, N. Abrosimov, A. Firsov, "A Graded Crystal Monochromator at BESSY II" SPIE 4145 , p.p. 122-128, (2000)

□

PERFORMANCE DATA

resolution: E/ΔE ~ 4 000	Exit flux: 10 ⁷ – 10 ¹⁰ phot/sec/100mA
spot size at experiment:	horizontal: 250 μm, vertical: 600 μm (5 μm x 5 μm with capillary optics)
Maximal divergence at sample position:	2.5 mrad (hor.) x 0.5 mrad (ver.)




TECHNICAL REFERENCE


FLOOR PLAN

in – air experimental arrangement

Top view




Side View



In – vacuum experimental arrangement.

Vacuum chamber position



Diagnostic Chamber: Hor./vert. slits 0-10 mm; vert. fixed slits 50 μm and 100μm; video monitor; filter assembly V, Ni, Fe, Cu foils, each of 1 μm thick; PIN diode.

GEOMETRICAL BOUNDARY CONDITIONS

In-air experimental area	Optical table 1,8m x 1.2m, in the hutch (see floor plan), 2 Capton windows, differential pumping.
Focus position in the hutch	Horizontal: 35 m from the source (see floor plan) Vertical: variable 30 m - 36 m from the source.
Focus height in the hutch	1465 - 1490 mm above floor (depending on lateral position)
In- vacuum chamber	User chamber must fit in 1.5 m overall length.
User flanges	CF35

Special instrumentation	<p>The beamline provides an experimental set-up for EXAFS and XANES measurements at-air and in He atmosphere in the energy range of 4 keV – 15 keV. Beam intensity is stabilized by MOSTAB electronics with an accuracy of 0.3 %. The detector system consists of three ionisation chambers, a Si-PIN photodiode for fluorescence measurements, a scintillation counter and an energy-dispersive detector (Röntec X-Flash). An add-on microprobe capillary system with the spatial resolution of $>5 \mu\text{m}$ is available for micro-EXAFS, micro-fluorescence and micro-diffraction experiments.</p> <p>For X-Ray diffractometry and reflectometry the experimental hutch is equipped with a 6-axis HUBER goniometer and a 2-dimensional cross-wire detector with a spatial resolution of $150 \mu\text{m}$. In addition detectors listed are available. A high-resolution X-ray CCD camera with a pixel size of $6.7 \mu\text{m}$ can be used for experiments alignment.</p>
-------------------------	---

□

VACUUM REQUIREMENTS

max. pressure	$\sim 6 \times 10^{-8}$ mbar at last valve for windowless experiment, live zero point signal for interlock
oil free vacuum system	yes
In - air experiments	yes

INFRASTRUCTURE AT EXPERIMENTAL STATION

electrical power supplies	220V, 380 V max 44 kVA
demineralized cooling water	closed system, forerun $30^\circ / 10$ bar, return run $\leq 42^\circ / 2$ bar
pressurized air	8 bar
oilfree exhaust line	not for hazardous gases!
He-recycling system	yes

DATA ACQUISITION

control system	PC-based BESSY monochromator control system EMP/2
data-acquisition computer	Personal Computer, measurement bus-extension, OS/2-operating system
data-acquisition software	Windows-NT, RADICON RDPW software, SPEC software for 6-axis HUBER diffractometer.
remote-control	V24-serial port, DMC/AMC-protocol