


POLARIMETER Soft X-Ray (UHV-) 8 Axes Diffractometer	 0049 030 8062 - 14900
--	---

CONTACT PERSON	Building	Room	Phone	e-mail
F. Schäfers	14.51	3313	12946	franz.schaefers@helmholtz-berlin.de
A. Sokolov	14.51	3370	12986	andrey.sokolov@helmholtz-berlin.de
M. Mertin	14.51	3354	12930	marcel.mertin@helmholtz-berlin.de



The high precision 8-axes ultra-high vacuum compatible (UHV)-polarimeter is a multipurpose instrument which can be used as a selfcalibrating polarisation detector for linearly and circularly polarised UV- and soft X-ray light [1, 2]. It can also be used for the characterisation of either reflection or transmission properties (reflectometer) or polarising and phase retarding properties (ellipsometer) of any optical element [3]. Magneto-optical experiments are possible in transmission, as the XMCD [4] or XMLD [5] (Magnetic Circular / Linear Dichroism) that are intensity measurements. Additionally a polarisation analysis of the transmitted light is possible which allows for Faraday- [4] or Voigt-measurements [5]. In reflection the magneto-optical Kerr effect can be exploited in longitudinal (L-MOKE) or transversal (T-MOKE) geometry as intensity measurement [6 - 8] to investigate thin films as well as magnetic multilayers [9]. Independent two-dimensional rotation of the detector enables any non-specular magnetic scattering experiment on magnetic dots or grains. A load-lock transfer chamber allows for quick and easy sample exchange.


Typical experiments

- **Reflectometry, ellipsometry (multilayer mirrors, films)**
 - Characterisation of optical elements
 - Reflection, transmission properties (s-, p-pol.)
 - Polarising properties (phase retardation)
- **Polarimetry (with two optical elements)**
 - Determination of polarisation of incident light (Stokes $S_{0,1,2,3}$)
 - Polarisers: transmission multilayers, $\lambda/4$ -plates
 - Linear Analysers: reflection multilayers, mirrors
- **Magneto-optical spectroscopy (in reflection, transmission)**
 - Resonant Magnetic Scattering (specular and diffuse)
 - Intensity spectroscopy: MCD, LMD, Kerr-effect (L, T-MOKE)
 - Polarisation spectroscopy: Faraday-, Voigt-effect


References

- [1] F. Schäfers et al., *A Soft Soft X-Ray Polarimeter using Multilayer Optics: Complete Analysis of the Polarization State of Light*, Appl. Opt. **38**, 4074-4088 (1999)
- [2] F. Schäfers, *Multilayers for the EUV/Soft X-ray Range*, Physica B **283**, 119-124 (2000)
- [3] H.-Ch. Mertins et al., *Observation of the X-ray Magneto-Optical Voigt Effect*, Phys. Rev. Lett. **87**, 047401-1-4, (2001)
- [4] H.-Ch. Mertins et al., *Resonant magnetic reflection coefficients at the Fe 2p edge obtained with linearly and circularly polarized soft x rays*, Phys. Rev. B **66**, 184404-1-8 (2002)
- [5] O. Zaharko et al., *Exchange coupling in Fe/NiO/Co film studied by soft x-ray resonant magnetic reflectivity*, Phys. Rev. B **66**, 134406 (2002)
- [6] Z. Wang et al., *Broadband multilayer polarizers for the extreme ultraviolet*, J. Appl. Phys. **99** (5) 056108-1-3 (2006)
- [7] F. Eriksson et al., *Interface Engineering of Short-period Ni/V Multilayer X-ray Mirrors*, Thin Solid Films **500**, 1-2, 84-95 (2006)
- [8] F. Schäfers, *Multilayer-based Soft X-ray Polarimetry*, Optics and Precision Engineering (2007)


Technical sketch of the UHV 8-axes configuration




EXAMPLES




Reflectivity curve of a multilayer mirror - in s- and p-polarisation geometry ((2:1; $\theta:2\theta$ angular scan) [3]




Determination of the complete polarisation state (Stokes vector) of incident light [1]




Faraday measurement on a thin Fe-film [4]



Phase retardation of a transmission multilayer [1]



T-MOKE, Reflectance $R_{T+/-}$ and corresponding asymmetry [6]



Hysteresis-curve in reflection

TECHNICAL DATA	
Io - light diagnostics	
Higher order filters	Be, B, C ₆ H ₈ , Ti, Cr, Fe
Collimator pinholes	Ø 0.2 – 2.0 mm
Detectors	Au-mesh, Fluorescence screen, GaAsP photodiode
Samples in transmission or reflection	polarisers, analysers
Maximum dimension	50 x 50 x 11 mm ³
Minimum dimension	10 x 10 x 0.5 mm ³
Incidence angle scan range	0° ≤ Θ _P , Θ _A ≤ 90°
Azimuthal angle scan range	0° ≤ α, β ≤ 370°
Minimum angle to normal incidence	4.5°
Minimum step size	0.001°

Heating	200° C
Magazine store	in-situ change of 10 samples
Load-lock transfer chamber	for 5 samples
Magnetic fields	$-450 \leq H \leq 450$ Oe
Trans.: in-/off-plane	
Refl: long./transv.	
Detector	GaAsP-photodiode with Keithley electrometer 617 (6514)
Dark current	$30 \cdot 10^{-15}$ A
Dynamic range	up to 8 orders of magnitude
Detector size	$4 \times 4 \text{ mm}^2$; $0.2 \times 4 \text{ mm}^2$
Scan range in plane	$0^\circ \leq 2\Theta_A \leq 180^\circ$
off-plane	$-10^\circ \Theta_D \leq +27^\circ$
Min. step size in plane	0.001°
off plane	0.001°
Sample – Detector Distance	150 mm
UHV-chamber	
Vacuum	$1 \cdot 10^{-8}$
Beam height	1400 ± 300 mm
Adjustments w.r.t. SR:	x, y, pitch, yaw, roll
Computer control	
Hardware (speckle.exp.bessy.de)	Intel(R) Pentium(R) D CPU
Software	3.40GHz, 1GB RAM, 250 GB Linux (Debian Etch) / SPEC 5.06
Scan options	h ν , θ , 2θ , detector off-plane, θ - 2θ azimuth α , β , polar- α , - β
multilayers, crystals:	Bragg-peak-scan (h ν - θ - 2θ)
gratings:	constant inc./diff./dev. angle